


CELL COUNTER

B CytoCube (small, compact, and accurate)

\$ P N Q B D U % F Ease of use / One-step Operation


- It takes only 2s to measure a sample.
- A slide can hold 6 samples.



Add sample



Insert the slide





One-click Analysis


Reliability

Designed for accurate and reliable cell counting

- 5-megapixel color CMOS for image capture
- Corresponding indicator turns on after sample introduction and ensures proper imaging and data capture.
- Optimized optical path provides excellent image quality.



Trypan Blue method to distinguish dead and living cells



Proven Cell Line

Excellent User Experience

User-friendly software interface

- Thorough detection of concentration, viability, diameter and roundness of total cells
- User level management system: Different administrative level permission maintains data confidentiality.
- Multiple formats for exporting data including PDF, EXCEL and JPG to meet different requirements of users.